



- ・どんな歯牙にもマッチする半透明色でできています。
- ・メカニカルロックベースの為、接着剤を問いません。

- ・丸みをおびたウィングは頬粘膜への刺激を軽減します。
- ・強度の高い多結晶アルミナ製です。

ロスタイプ	Torque	Ang.	Right/Left	.018 スロット	.022 スロット
  <u>1   1</u>	+12°	+5°	R L	<b>18-951R</b> <b>18-951L</b>	<b>22-951R</b> <b>22-951L</b>
  <u>2   2</u>	+8°	+9°	R L	<b>18-952R</b> <b>18-952L</b>	<b>22-952R</b> <b>22-952L</b>
  <u>3   3</u>	-2°	10°	R L	<b>18-953R</b> <b>18-953L</b>	<b>22-953R</b> <b>22-953L</b>
  <u>54   45</u>	-7°	0°	R L	<b>18-954R</b> <b>18-954L</b>	<b>22-954R</b> <b>22-954L</b>
  <u>21   12</u>	-1°	0°	R L	<b>18-961</b>	<b>22-961</b>
  <u>3   3</u>	-11°	+7°	R L	<b>18-963R</b> <b>18-963L</b>	<b>22-963R</b> <b>22-963L</b>
  <u>4   4</u>	-17°	0°	R L	<b>18-964R</b> <b>18-964L</b>	<b>22-964R</b> <b>22-964L</b>
  <u>5   5</u>	-22°	0°	R L	<b>18-965R</b> <b>18-965L</b>	<b>22-965R</b> <b>22-965L</b>
					1pk 10個入
1 症例キット	<u>5~5</u> <u>5~5</u>			<b>18-950K</b>	<b>22-950K</b>  1kit 20個入